


SEIZURE DETECTION USING A WEARABLE SENSOR

Team FitByters

Introduction

What is a Seizure?

- A sudden, uncontrolled electrical disturbance in the brain.
- Disruption in Electrical Pulse Patterns of the brain.
- Brain cells moving simultaneously at a faster rate than normal.
- Spreading to different parts of the brain, leading to sensations, awareness and movement changes.
- Maybe severe or subtle.
- May lead to SUDEP



Project Objectives

Detecting and Avoiding Seizures and Preventing SUDEP

- Design a system that can detect the event of seizures in an epilepsy patient in real-time.
- 2. Relay the information via an application to dear ones.
- Transmit the GPS location of the epilepsy patient in the event of a seizure.


System Devices



FitByte App







The "Brain"



4G LTE Base shield HAT



ASUS USB-BT400 w/Bluetooth Dongle




Antenna - 100mm



WitMotion BWT901CL Bluetooth mpu9250 + Gyro + Accelerometer(Triple Axis) + Angle + Magnetometer





Quectel EC25 Mini PCIe 4G/LTE Module




Raspberry Pi 4


FitByte





GPS Location






System Architecture

The Seizure Detection System comprises:


- 1. Bluetooth Gyro Sensor
- 2. Asus Bluetooth Transceiver
- 3. Raspberry Pi 4 Processor
- 4. 4G/LTE HAT:
 - a. LTE Base HAT
 - b. QUECTEL EC25 Mini PCIe 4G/LTE Module
 - c. Antenna
- 5. FitByte





Acceleration samples processing technique


```
Input acceleration
              samples
If (current_sample>previous_sample){
output value=+1
Elseif (current_sampleprevious_sample){
                                                     A "SWITCH" is defined as a 1 to -1
output value=-1
                                                     or a -1 to 1 transition.
Else{
output value=0
```


Switch samples













Finite State Machine of the System

Flowchart





Seizure



True Negative case



FitByte



System Demonstration

Acknowledgement

- 1. Dr. Yuhua Chen
- 2. Shreyas Poyrekar

References

- 1. Lliwen A Jones and Rhys H Thomas. Sudden death in epilepsy: Insights from the last 25 years. Seizure, 44:232–236, 2017.
- MayoClinic.Seizures.https://www.mayoclinic.org/diseases-conditions/seizure/symptoms-causes/syc-20365711, 2019. Online; accessed 12 November 2019.
- 3. Yash Paul. Various epileptic seizure detection techniques using biomedical signals: a review.Brain informatics, 5(2):6, 2018.

Thank You

